Med-Fit

Wireless 4 Channel TENS Machine

CE 2862

Patient Instructions & User Manual

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email: sales@med-fit.co.uk

www.med-fit.store

GENERAL DESCRIPTION	3
CONTRAINDICATIONS.....	4
PRECAUTIONS	4
WARNINGS.....	5
INTRODUCTION TO TENS	7
CONTENTS & GENERAL INFORMATION	8
OPERATING INSTRUCTIONS	9
REMOTE CONTROL FEATURES	11
MODULE FEATURES.....	12
CHARGING THE REMOTE CONTROL	13
CHARGING THE MODULES.....	14
INSTRUCTIONS FOR USE.....	15
ELECTRODE INSTRUCTIONS FOR USE	17
TROUBLESHOOTING	19
MAINTENANCE, TRANSPORTATION AND STORAGE	20
HELPFUL TIPS.....	21
CONFORMITY TO SAFETY STANDARDS	22
TECHNICAL SPECIFICATIONS	22
GLOSSARY OF SYMBOLS	23
ELECTROMAGNETIC COMPATIBILITY DESCRIPTIONS.....	

GENERAL DESCRIPTION

The Four Channel TENS is a battery operated single channel electronic nerve stimulator.

The wireless TENS is a four channel transcutaneous electrical nerve stimulator (TENS). The device generates small pulses of electrical current and delivers the pulse to the user through self-adhesive electrodes such that the underlying nerves are activated, and the pain associated with sore and aching muscles is relieved. The Four Channel TENS is portable, compact, rechargeable and stylish in design. It may help relieve pain in various areas of the body, such as lower back, shoulders, arms and legs.

The device generates small pulses of electrical current and delivers the pulse to the user's skin through adhesive gel electrode pads such that the underlying nerves are activated and the pain associated with sore and aching muscles is relieved

The TENS is portable, compact, rechargeable and stylish in design. It may help to relieve pain in various areas of the body, such as the lower back, shoulders, arms and legs.

INTENDED USE

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, arm, leg and foot, due to strain from exercise or normal household and work activities.

CONTRAINDICATIONS

1. Patients with pacemakers and cardiovascular problems should not use the device.
2. The device should not be administered over the abdomen or pelvis during pregnancy.
3. Patients with epilepsy should not apply the electrodes to the neck or head.
4. Patients with heart disease, cancer or any other health conditions should not use the device, unless recommended by a medical practitioner.
5. Patients who have suffered acute trauma or surgical procedures in the past six months.
6. Patients with undiagnosed pain syndromes should not use the device, unless recommended by a medical practitioner.

PRECAUTIONS

1. Do not apply the electrodes on inappropriate sites, such as over chest, the head, anterior neck, eyes, mouth and the front of throat.
2. Do not apply the electrodes on areas of broken or damaged skin, such as open wounds.
3. Do not apply the electrodes on the swollen, infected or inflamed areas or skin eruptions.
4. Do not apply the electrodes on the skin with diminished sensation because nerve damage is likely to diminish the effectiveness and the patient maybe unaware that high-intensity currents are causing skin irritation.
5. Patients should not use the device in water or when operating hazardous machinery such as driving.
6. Do not use the device near monitoring equipment.
7. Do not use the device close to transdermal drug delivery system.

WARNINGS

1. PATIENT with an implanted electronic device (such as a cardiac pacemaker) should not use the device unless specialist medical opinion has first been obtained.
2. Simultaneous connection of a PATIENT to a high frequency surgical ME EQUIPMENT may result in burns at the site of the STIMULATOR electrodes and possible damage to the STIMULATOR.
3. Operation in close proximity (e.g. 1m) to a shortwave or microwave therapy ME EQUIPMENT may produce instability in the STIMULATOR output.
4. Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
5. Stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.
6. Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
7. Stimulation should not be applied transthoracically in that the introduction of electrical current into the heart may cause cardiac arrhythmias.
8. Stimulation should not be applied transcerebrally.
9. Never use the device while driving, operating machinery or during activities in which involuntary muscle contractions may endanger the user or others.
10. Do not use the device in the bath or shower.
11. Never use the device while sleeping.
12. Never immerse the device in any liquid.
13. No service and no maintenance shall be done while in use with a patient.

WARNINGS

14. Do not modify the device without authorization of the manufacturer.
15. Keep unit out of the reach of young children.
16. Apply the whole surface of the electrodes firmly to the skin. Do not use electrodes that do not stick properly to the skin or only partially stick to the skin.
17. The electrodes should be discarded when they are no longer adhering.
18. Only use the electrodes provided by the manufacturer, do not use other kinds of electrodes. Otherwise, it will cause skin irritation if the current density is more than $2\text{mA}/\text{cm}^2$.
19. Over long application may cause harm to the patient.
20. Stimulation should not be done in conjunction with OXYGEN RICH ENVIRONMENTS.

ADVERSE REACTIONS

1. You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin;
2. You should stop using the device and should consult with your healthcare advisor if you experience adverse reactions from the unit.

INTRODUCTION TO TENS

What is TENS?

Transcutaneous electrical nerve stimulation is a pain control treatment. It is often called TENS for short.

A TENS unit is a portable, pocket-sized, battery-powered device.


The TENS unit uses mild, safe electrical signals to help control pain and delivers the electrical signal to the body through self-adhesive conductive electrodes.

How does TENS work?

The most common TENS programmes use high-frequency stimulation, which is the first choice for both acute and chronic pain. High-frequency stimulation sends impulses to the nervous system's own pain-inhibiting mechanisms, which block the pain. You can use it as often and as long as you like, but each treatment should last at least 1 hour.


Another type of TENS is low-frequency stimulation. Low-frequency TENS treatment can alleviate pain by stimulating muscles to release the body's own morphine-like substances, endorphins.

Place the electrodes on a muscle in the painful area so that a visible contraction occurs.

CONTENTS & GENERAL INFORMATION

Please check carefully the contents of the Four Channel Wireless TENS and Accessories

1. Painless TENS modules x4

2. Remote control x1

3. Large lumbar electrodes x2

4. Self adhesive electrodes x8

5. Butterfly electrodes x2

6. AC adaptor x1


7. Quad charging lead x1

8. Instruction & user manual x1

Re-ordering self adhesive electrodes

RC-400 20cm X 10cm RC-360 5CM Round RC-240 15cm x 9cm

OPERATING INSTRUCTIONS


No Function description

1. Programme Display
2. Intensity Display
3. Intensity Up
4. Intensity Down
5. On/Off Button and Mode Selector
6. Programme Selector
7. Module Channel Display Screen
8. Treatment Timer
9. Battery Indicator

10 OPERATING INSTRUCTIONS

The TENS Module connected to the butterfly electrode

Module connected to the round electrode pad.

FIG B

When the TENS is switched on and not in use, it will automatically turn off after approximately one minute, to conserve battery life.

Please note when using the small round shaped electrodes (see Fig B) you must connect 2 electrodes, one to each connection point.

OPERATING INSTRUCTIONS

Remote control features


On/Off Button & Module selection. [Button 5].

 Press the On/Off button once to turn on the remote control, now press the On/Off button again to select the required module 1, 2, 3, or 4. This is displayed in the bottom right segment of the remote control display. To turn off your remote control, press and hold the On/Off button for two seconds.

Programme selector button

Each single press of the [P] button changes the setting to the next programme you may set any programme for each module. [Button 6].

Programme 1
Press P to change the programmes


90 minutes timer
Module 1
Press the On/Off button to change module selection


Detailed description: A diagram of a remote control display. The display is divided into four segments. The top-left segment shows a 'P' button icon. The top-right segment shows '1' and '90'. The bottom-left segment shows '00' and a signal strength indicator. The bottom-right segment shows '1' and a clock icon. Green arrows point from the text labels to the corresponding segments: 'Programme 1' to the '1', '90 minutes timer' to the '90', 'Module 1' to the '1', and 'Press the On/Off button to change module selection' to the clock icon.

For example: You can use the same programme on each module or you can set a different programme for each module as shown below.

Programme 1
Module 1

Programme 4
Module 2

Examples of possible settings

Detailed description: Two diagrams of remote control displays. The left display shows '1' and '90' in the top-right, '00' in the bottom-left, and '1' in the bottom-right. The right display shows '4' and '90' in the top-right, '00' in the bottom-left, and '2' in the bottom-right. Green arrows point from the text labels to the corresponding segments: 'Programme 1' to the '1', 'Module 1' to the '1', 'Programme 4' to the '4', and 'Module 2' to the '2'.

12 OPERATING INSTRUCTIONS

Intensity settings


This button increases and decreases the intensity for each module which allows you to set individual levels for each module ideal when treating multiple areas of the body [Button 3 and 4].

The intensity is adjustable from level 1 to level 16 for each module.


Module features

The four modules are marked as shown.

A green light indicates the module is turned on.

On/Off button

Press stud female connectors

Note: When using the individual 5cm round electrodes, both female connectors must be connected to the electrode pads. You must connect two pads to the module.

To turn on the module, press the On/Off button once a single beep and the green light indicate the module is on. A further press of the On/Off button turns the module off.

CHARGING THE REMOTE CONTROL

13

It is always recommended that you charge your remote control before using it for the first time.

Charging indicator flashes when charging.

There are two methods of charging your TENS remote control.

To charge from a mains outlet socket, simply connect the USB cable to the mains charging adaptor and the other end to the remote control. Please do not force the USB cable into the sockets see Fig C

Charging the TENS Remote via a computer USB port.

You may charge your TENS remote control through a computer by using the USB cable port see Fig D

Fig C

Fig D

USB Charging Port

CHARGING THE MODULES

It is always recommended that you charge your TENS modules prior to using it for the first time.

There are two methods of charging your TENS Modules.

To charge from a mains outlet socket, simply connect the USB cable to the mains charging adaptor and the other end to the TENS Module. Please do not force the USB cable into the sockets see Fig C

Charging the TENS Module via a computer USB port.

You may charge your TENS Module through a computer by using the USB cable port see Fig D

How long should I charge the TENS Module.

Once connected the charging light will flash green a fully charged battery is indicated by a continuous green LED.

Auto shut down mode.

To save battery life the TENS will automatically shut down once switched on but is not being used.


Fig C

Fig D

INSTRUCTIONS FOR USE


15

Step by Step

For your very first treatment session, it is recommended that you use just one module so that you familiarise yourself with the operation.

Please follow the steps below.


1. Make sure both inputs of the TENS Module are connected into the snap-pin connectors of the gel electrode pad or pads.
2. Peel off the plastic film from the opposite side of the butterfly or standard gel electrodes supplied.
3. Place the TENS module attached to the gel electrode pad onto the treatment area. Press down firmly and ensure a full and firm contact is made with the skin.
4. Press the On/Off power button to turn on the TENS module as indicated by the LED light being green.
5. Turn on your remote control you should hear two beeps this indicates the remote and tens modules are responding. We recommend starting with P1 as this is ideal for first-time use.
6. The remote default setting is P1 and Module 1. Please use these settings for your first treatment.


7. Increase the intensity by pressing the [+] button each press increases the intensity and is indicated by a single audible beep. Once you increase the intensity, it will count down in minute intervals. Set a level of stimulation that is a pleasant sensation and not too strong. On this setting you should feel a constant tingling sensation.

Important notes to remember


Always turn on the module or modules before turning on the remote control. Always check the remote control channel is correct for the module you are using.

Make sure the number on the remote matches the module.

If the number on the remote does not match the number on the module it will not change the setting, make sure the numbers match.

ELECTRODE INSTRUCTIONS FOR USE

Turn Stimulator OFF before applying or removing electrodes

Application of the self-adhesive pads & aftercare

The electrode pads have a maximum usage of 25-30 applications to keep your pads in good condition apply a small amount of water to the gel side of the pad every 4th application (do not over wet). Do not apply any gels or creams to the pads as this will destroy the hydro-gel on the electrode pad.

Application

1. Skin site must be very clean and dry, dirty, flaky or oily skin will prevent electrodes from adhering to the skin. If necessary, trim excess hair with scissors. If skin is oily wipe down with an alcohol or electrode skin prep prior to application. Be sure to wash hands before handling electrodes.
2. Remove electrodes from bag and reseal bag to protect remaining electrodes.
3. Insert the stud electrode on to the back of the device.
4. Grasping a tiny edge of the electrode, peel and remove electrode from the protective plastic liner. Save liner for electrode storage.
5. Place electrode onto skin treatment site (as recommended by your clinician) by firmly applying from the center of the electrode to the outer edges. Adhesion improves when electrodes reach skin temperature.
6. If gel appears over-saturated with excessive moisture or perspiration, allow the electrode to air dry in a refrigerator with the gel side facing up until the gel regains its tack. If the gel appears dry, try adding a few drops of water to the gel and allow to rest in a dust-free environment until the gel regains its tack.

It is important to smear water on to the gel of the electrode every 3rd or 4th application as this increases the adhesion of the gel.

Removal and storage

1. Lift a corner of the electrode and slowly peel the electrode off the skin, touching the adhesive gel as little as possible.
2. Place the electrodes back onto the saved protective plastic liner.
3. Return the electrodes back into the storage bag and reseal tightly to prevent dry-out.
4. Store at room or cool temperature and keep out of direct sunlight.
5. The life of the electrode varies depending on skin conditions, amount of use, storage and climate. Electrode life may be extended by carefully following the application, removal, and storage instructions.

Caution

1. DO NOT place electrodes on broken skin. If skin irritation develops discontinue use. Consult physician. Replace electrodes when they do not adhere or when treatment becomes uncomfortable.
2. DO NOT use unit while driving or operating machinery
3. DO NOT wear electrodes when showering, bathing or swimming
4. DO NOT apply electrodes across the head or across the heart or on the front neck.
5. Keep electrodes separated during treatment
6. DO NOT exceed 0.1 watts/cm²
7. Using stimulation electrodes that are small or incorrectly applied could result in discomfort or skin burns.

Latex Free

When the TENS is switched on, pushing the [+] button, if the device emits 1 bleep and the green light flashes once, but you feel no sensation through the electrode pads this may be due to the following conditions:

1. The TENS Pads are not securely fastened to the skin.
 2. The area being treated is too dry please remove your pads and rub warm water into the treatment area and reapply the pads.
 3. The electrode pads need replacing.
- * Please note self-adhesive pads need replacing after 20-30 applications. To increase the life of your pads smear small amount of water onto the pad every 4th application.

If your device is not operating properly please check for common problems and suggested solutions. If the recommended action does not solve the problem, please contact the seller.

Stimulation is weak or non-existent

- Be sure skin is clean and gel electrode pads are firmly attached to the skin.
- The battery is low and needs to be charged.
- The gel electrode pad needs to be replaced.

1. Turn off the device before cleaning.
2. The electrode pad of the device can be cleaned with a soft and wet cloth.
3. Never use excessive water or cleaning agents to clean the device.
4. The gel electrode pad can remain on the control unit between uses, and are covered by the plastic film to avoid getting dirty and dry.
5. Avoid exposure to extreme moisture or temperature.



HELPFUL TIPS

Once you have familiarised yourself with the controls and features of your TENS device, it is important to place the TENS electrodes in a position which gives the most pain relief. This may take 3 or 4 attempts to find the most suitable position, for maximum pain relief.

If you are using two electrodes, place the electrodes directly onto the painful area at a position where you feel the pain starts and where it finishes. You may now position the electrodes around the painful area to locate the most suitable position for maximum pain relief.

The alternative method is to use four electrodes surrounding the painful area see examples below

The complete area between the electrodes will now be treated when positioning the electrodes as shown.

Statement of EMC

IEC 60601-1-2:2014, EN 60601-1-2:2015, IEC 60601-1-11:2015 clause 12, EN 60601-1-11:2015 clause 12, IEC 60601-2-10:2012 + A1:2016 clause 201.17 & 202, EN 60601-2-10:2015 + A1:2016 clause 201.17 & 202

Conformity to MDD Requirements

IEC 60601-1:2005+A1:2012, EN 60601-1:2006+A1:2013, ANSI/AAMI ES60601-1:2005+A1:2012, C1:2009 and A2:2010, IEC 60601-1-11:2010, EN 60601-1-11:2010, IEC 60601-2-10:2012, EN 60601-2-10:2015

TECHNICAL SPECIFICATIONS

Transcutaneous Electrical Nerve Stimulators: KTR-2492

Operating condition Temperature: +5°C+40°C

Relative humidity: 15% - 93%RH

Atmospheric pressure: 700hPa-1060hPa

Environment for storage

Temperature: -25°C+70°C

Humidity: 0-93%RH

Barometric pressure: 700hPa - 1060hPa

Shenzhen Kentro Medical Electronics Co.,Ltd
2nd Floor, No. 11, Shanzhuang Road, Xikeng Village,
Yuanshan Street, Longgang District, Shenzhen City,
Guangdong Province, China


WellKang Ltd (www.CE-marking.eu) Enterprise Hub,
NW Business Complex, 1 Beraghmore Road, Derry,
BT48 8SE, Northern Ireland.

WellKang Ltd (www.CE-marking.eu) Enterprise Hub,
NW Business Complex, 1 Beraghmore Road, Derry,
BT48 8SE, Northern Ireland.

GLOSSARY OF SYMBOLS

	Authorized Representative in the UK
	Authorized Representative in the European Community
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Class II equipment
	Date of manufacture
	Manufacturer
	Specifies serial number
	Type BF applied part
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Follow instructions for use.
IP22	The first number 2: Protected against solid foreign objects of 12.5 mm Φ and greater. The second number: Protected against vertically falling water drops when enclosure is tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.

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ELECTROMAGNETIC COMPATIBILITY DESCRIPTIONS

This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

- 1) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 2) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 3) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emission		
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of DEVICE should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The DEVICE is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations flicker emissions IEC 61000-3-3	Complies	

1.1 Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input /output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ± 2 kV common mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % U_T ; 250/300 cycle	0 % U_T ; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % U_T ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEVICE requires continued operation during power mains interruptions, it is recommended that the DEVICE be powered from an uninterruptible power supply or a battery.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a. c. mains voltage prior to application of the test level.			

1.2 Recommended isolation distance

Guidance and manufacturer's declaration – electromagnetic immunity –for EQUIPMENT and SYSTEM

Guidance and manufacturer's declaration – electromagnetic immunity			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz	3V 150 kHz to 80 MHz	

	<p>6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz</p> <p>10 V/m</p> <p>80 MHz to 2.7 GHz</p> <p>385MHz-5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)</p>	<p>6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz</p> <p>10 V/m</p> <p>80 MHz to 2.7 GHz</p> <p>385MHz-5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Model + Name, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>1.2.1 Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P} \quad d = \left[\frac{12}{V_2} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.7 \text{ GHz}$ <p>where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	---	---	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

- a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- b Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEVICE is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DEVICE.
- c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the DEVICE

The DEVICE is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the DEVICE can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DEVICE as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter W	Separation distance according to frequency of transmitter			
	50 kHz to 80 MHz outside ISM and amateur radio bands	150 kHz to 80 MHz in ISM and amateur radio bands	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{12}{V_2}\right]\sqrt{P}$	$d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	$d = \left[\frac{7}{E_1}\right]\sqrt{P}$

0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.035	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	35	70

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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EC REP

WellKang Ltd (www.CE-marking.eu) Enterprise Hub,
NW Business Complex, 1 Beraghmore Road,Derry, BT48 8SE, Northern Ireland.

UKRP

WellKang Ltd (www.CE-marking.eu) Enterprise Hub,
NW Business Complex, 1 Beraghmore Road,Derry, BT48 8SE, Northern Ireland.

Med-Fit UK Ltd.

Unit 8 Martel Court, S. Park Business Park, Hamilton Road, Stockport, SK1 2AF

Tel: 0161 429 7330 | Tel: 0161 429 7330 | Fax: 0161 427 0215

Email: sales@med-fit.co.uk | www.tensmachineuk.com

Company registration number 08758741 | Vat registration number 308286105